

MASTERWORKS™ Ortho – Instructions for Use



Master Works™ Ortho is supplied in a duo pot containing the two components A and B.



The duo pot contains sufficient material, when mixed together, for two castings, a total quantity of 1066cc.



Firstly remove inner top pot from the larger main pot.



Remove protective film from the A component pot.



Remove lid from inner, B component pot.



Remove protective film from the B Component.



Pour contents of the B component pot into main A component pot. Ensure all the B component is used.



Mix the two components together with a spatula until a uniform color is achieved.



Alternatively use a Master Works™ small high shear mixing blade in a hand held electric drill.

MASTERWORKS™ Ortho – Instructions for Use (cont.)



Mix the two components together until a uniform color is achieved.



Remove the mixing blade from the pot.



Wash the mixing blade by spinning in water.



Pour the mixed material from a height of about 300 mm into one end of the impression.



Allow the material to flow forward to fill the impression.



Fill the second impression as before.



Give a gentle shake to the impression box to ensure material is level.



Allow the filled impressions to set. Initial set if a few minutes. De-mold circa 30 minutes.



Following cure remove any surplus foam with a dry, stiff brush.

MASTERWORKS™ is a licensed trademark. These suggestions and data are based on information we believe to be reliable. They are offered in good faith, but without guarantee, as conditions and methods of use of our products are beyond our control. We recommend that the prospective user determine the suitability of our material and suggestions before adopting them on a commercial scale. Suggestions for uses of our products or the inclusion of descriptive material from patents and the citation of specific patents in this publication should not be understood as recommending the use of our products in violation of any patent or as permission or license to use any patents of PUMA Polymers. March 25, 2009